

WELCOME

TO



FDSPharmacy

Learn and Educate

This is an Education Platform

We provide Free PDF Notes and Videos Classes for Pharmacy Students

Web Site <http://www.fdspharmacy.in/>

You tube <https://www.youtube.com/channel/UC77iEsiuZolU4pB8WAJIR5Q>

What app <https://chat.whatsapp.com/IzSgXtFEvhS4LN5xhUgq5z>

Telegram <https://t.me/+cvxm17xSloA4MjVI>

Face book <https://www.facebook.com/Fdspharmacy-105764311994440/>

E-mail fdspharmacyinfo@gmail.com

Diploma in Pharmacy 1st Year

Pharmaceutics

Experiment

To demonstrate the friability test of the tablets as per the monographs.

Aim:

To demonstrate the friability test of the tablets as per the monographs.

Reference :

‘ Dr. Gupta G.D , Dr. Sharma Shailish , Dr. Sharma Neelam ’
“Practical Manual of Pharmaceutics” Published by Nirali Prakashan, Page
no 136 – 139

Apparatus and Materials Required :

Tablets and Roche friabilator.

Theory :

Friability test is performed to evaluate the ability of the tablets to withstand abrasion in packing handling and transporting. The instrument used for this test is known as friability test apparatus or friabilator.

Friability of a tablet can be determined in laboratory by Roche Friabilator (figure 12). This consists of a plastic chamber that revolves at 25rpm. The tablets are dropped from a distance of 6 inches in the friabilator which is then operated for 100 revolutions. The tablets are reweighed and friability (F) is calculated as:

$$F = \frac{W_{\text{initial}} - W_{\text{final}}}{W_{\text{initial}}}$$

Tablet that loses less than 0.1-0.5 % of the tablet weight is considered acceptable.




Figure 12: Roche Friabilator


Procedure

1. The cleanliness of the equipment should be checked.
2. The instrument should be switched on, drum should be initialised itself to the loading position, the display should now show START.
3. Weight of the tablet should be weighed and recorded..
4. 20 tablets should be taken if individual tablet weight is up to 0.650 g and 10 tablets should be taken if individual tablet weight is above 0.650 g.
5. The counts should be adjusted to 100 by pressing the COUNT key followed by '1', '0', '0' keys (Count range 1 to 99999).
6. ENTER key should be pressed to confirm the reading.
7. COUNT key should be pressed to see and confirm the number of counts.
8. On confirming the number of counts RUN/HALT key should be pressed to start, the display that shows the elapsed count.
9. Test can be performed by adjusting "TIME similarly as "COUNT".
10. When test is over drum should rotate in reverse direction for discharging the tablets in the tray.
11. The test over should be indicated by an audible beep and display should show END.

12. The drum should initialise itself to loading position and display should show START indicating the instrument is ready for the next run.
13. The tablets should be removed from the tray.
14. The tablets should be de-dusted and weighed and the weight of the tablets should be noted down.
15. The percentage loss in the weight should be calculated by using the formula.
16. The drums should not be held while rotating.
17. The tablets should be destroyed after testing.

Result :

The friability test of the tablets as per the monographs was performed.



Hello

Friends

**If you Get Any Help From This Notes /
Videos**

Next You Turn To Help Us

Please Contribute Some Amount

To Our

FDSPharmacy Team

Phone pe 6398439940

Paytm 6398439940

Google Pay 6398439940



Amir Khan

Thank You

Keep Supporting