

# WELCOME

# TO



# FDSPharmacy

Learn and Educate

This is an Education Platform

We provide Free PDF Notes and Videos Classes for Pharmacy Students

Web Site <http://www.fdspharmacy.in/>

You tube <https://www.youtube.com/channel/UC77iEsiuZolU4pB8WAJIR5Q>

What app <https://chat.whatsapp.com/IzSgXtFEvhS4LN5xhUgq5z>

Telegram <https://t.me/+cvxm17xSloA4MjVI>

Face book <https://www.facebook.com/Fdspharmacy-105764311994440/>

E-mail [fdspharmacyinfo@gmail.com](mailto:fdspharmacyinfo@gmail.com)

# Diploma in Pharmacy 2<sup>nd</sup> Year

## Biochemistry & Clinical Pathology

### Experiment

To determine the urea in blood /serum.

#### Aim:

To determine the urea in blood /serum.

#### Reference :


‘ Dr. Gupta G.D. , Dr. Sharma Shailesh, Kaur Manpreet, “Practical Manual of Biochemistry & Clinical Pathology” Published by Nirali Prakashan, Page no 53 – 58

#### Materials Required

Distilled water, blood serum, photometer, pipette and cuvette.

#### Theory :

Liver converts ammonia (primarily produced by the breakdown of amino acids) into urea. These reactions are used in kinetic enzymatic estimation of urea:



Urea is hydrolysed by urease into ammonia. Ammonia and 2-oxoglutarate are combined by glutamate dehydrogenase to produce glutamate. A decrease in absorbance at 340 nm (Warburg's optical test) is used to photometrically detect the conversion from NADH to NAD in this reaction.

## Procedure

- 1) The photometer should be switched on and allowed to heat for 10 minutes at  $37^{\circ}\text{C}$ .
- 2) The wavelength should be set at 340 nm and distilled water should be used to make the blanking. All of the absorbance that is further described is read using distilled water as a reference.
- 3) 3 cuvettes are available, i.e., one for blank reaction, one for the standard reaction and one for the sample (blood serum) reaction.
- 4) **For Blank Reaction:** 0.02ml of distilled water and 2ml of the reagent (working solution) should be pipetted out into the cuvette. The stop-watch should be pressed at this moment. The solution should be pipetted out inside the cuvette once again up and down to ensure proper mixing. The cuvette should be transferred immediately into the heated photometer. The initial absorbance  $A_1$  should be measured exactly 30 seconds after the working solution has been pipetted, and the second absorbance  $A_2$  should be measured after 1 minute.

$A_1$ blank	
$A_2$ blank	

- 5) 0.02ml of distilled water and 2ml of the reagent (working solution) should be pipetted out into the cuvette. The stop-watch should be pressed at this moment. The solution should be pipetted out inside the cuvette once again up and down to ensure proper mixing. The cuvette should be transferred immediately into the heated photometer. The initial absorbance  $A_1$  should be measured exactly 30 seconds after the working solution has been pipetted, and the second absorbance  $A_2$  should be measured after 1 minute.

$A_1$ Standard	
$A_2$ Standard	

6) 0.02ml of the serum sample and 2ml of reagent (working solution) should be pipetted out.  $A_1$  and  $A_2$  should be measured as compared to the previous measurements.

$A_1$ Sample	
$A_2$ Sample	

7) The differences in absorbance should be calculated.

$\Delta A$ blank = $A_1$ blank - $A_2$ blank	
$\Delta A$ Standard = $A_1$ Standard - $A_2$ Standard	
$\Delta A$ Sample = $A_1$ Sample - $A_2$ Sample	

**Calculate the Concentration:**

Urea (mmol/l) =  $\Delta A$  Sample -  $\Delta A$  blank /  $\Delta A$  Standard -  $\Delta A$  blank. C standard (15 mmol/l)

**Result :**

Urea concentration in blood /serum was determined.

**Hello**

**Friends**

**If you Get Any Help From This Notes /  
Videos**

**Next You Turn To Help Us**

**Please Contribute Some Amount**

**To Our**

**FDSPharmacy Team**

**Phone Pe            6398439940**

**Paytm                6398439940**

**Google Pay        6398439940**



Amir Khan

**Thank You**

**Keep Supporting**