

This is an Education Platform

We provide Free PDF Notes and Videos Classes for Pharmacy Students

Web Site http://www.fdspharmacy.in/

You tube https://www.youtube.com/channel/UC77iEsiuZolU4pB8WAJIR5Q

What app https://chat.whatsapp.com/IzSgXtFEvhS4LN5xhUgq5z

Telegram https://t.me/+cvxmi7xSloA4MjVl

Face book <u>https://www.facebook.com/Fdspharmacy-105764311994440/</u>

E-mail fdspharmacyinfo@gmail.com



Diploma in Pharmacy 1st Year Pharmaceutical Chemistry Experiment

To perform the identification test and test for purity of caffeine. Aim:

To perform the identification test and test for purity of caffeine.

Reference :

⁶ Dr. Gupta G.D., Dr. Sharma Shailish, Kaur Baljeet ² "Practical Manual of Pharmaceutical Chemistry" Published by Nirali Prakashan, Page no 73 - 75

Requirements:

Test tube, porcelain dish, water bath, caffeine, conc. HCI, potassium chloride. bromomethylene blue, tannic acid solution, ammonia solution, and bromomethylene blue solution.

Procedure: Identification Tests

| S.No. | Experiments | Observations | Inferences |
|-------|---------------------------------------|-------------------|-------------|
| 1 | Add few drops of tannic acid | White precipitate | Presence of |
| | solution to a saturated solution. | is formed | Caffeine |
| | | | |
| 2 | In a porcelain dish, pure crystals of | A purple colour | Presence of |
| | caffeine are combined with | is formed. | Caffeine |
| | potassium chloride and conc. Add | | |
| | 2-3 drops of HCI and evaporated to | | |
| | dryness using a water bath and | | |
| | ammonia solution. | | |



Test for Purity

- Test for Heavy Metals: Ig of Caffeine should be dissolved in 5ml of 0.1ml HCI and diluted to 25ml with water. The solution should pass the heavy metals method A limit test (20ppm)
- 2. Test for Sulphated Ash: Not more than 0.1% determined by the method in 1.P.
- 3. Acidity/Alkalinity: 0.2g should be dissolved in 10ml boiling water and then cooled. 0.1ml of bromomethylene blue solution should be added. The solution is green/yellow coloured. 0.02M NaOH should be titrated to obtain blue colour, not more than 0.1ml is required.

Result: The identification test and test for purity of caffeine were performed as per I.P.



Hello

Friends

If you Get Any Help From This Notes /

Videos

Next You Turn To Help Us

Please Contribute Some Amount

To Our

FDSPharmacy Team

Phone pe 6398439940

Paytm 6398439940

Google Pay 6398439940



Amir Khan



