WELCOME



This is an Education Platform

We provide Free PDF Notes and Videos Classes for Pharmacy Students

Web Site http://www.fdspharmacy.in/

You tube https://www.youtube.com/channel/UC77iEsiuZolU4pB8WAJIR5Q

What app https://chat.whatsapp.com/IzSgXtFEvhS4LN5xhUgq5z

Telegram https://t.me/+cvxm17xSloA4MjVl

Face book https://www.facebook.com/Fdspharmacy-105764311994440/

E-mail fdspharmacyinfo@gmail.com



Diploma in Pharmacy 1st Year Human Anatomy & Physiology Experiment

Recording of pulse oxygen (before and after exertion).

Aim:

Recording of pulse oxygen (before and after exertion).

Reference:

Dr. Gupta G.D , Dr. Sharma Shailesh , Dr. Sharma Rahul Kumar , "Practical Manual of Human Anatomy and Physiology" Published by Nirali Prakashan , Pg.No 106 - 109

Apparatus Required:

Pulse oximeter

Theory:

A pulse oximeter is a device that measures the amount of oxygen in the blood (oxygen saturation). A pulse oximeter can be used to calculate blood oxygen levels. It calculates the blood oxygen saturation and pulse rate using laser beams. The amount of oxygen carried in the blood is determined by oxygen saturation Without taking a blood sample, the pulse oximeter may determine the amount of oxygen in the blood.

Procedure

- 1) Any nail polish/false nails should be removed and hand should be warned (if cold).
- **2)** Before taking measurement, resting should be done for at least 5 minute.
- 3) The hand should be relaxed at the chest at heart level and held still
- **4)** Oximeter should be switched on and placed on the middle or index finger.
- 5) The oximeter should be kept in place for at least one minute or longer if the reading is not stable as the reading may take time to steady.

- **6)** The highest result once it has not changed for 5 seconds should be noted.
- 7) Each reading should be identified carefully.
- **8)** Recording should be started from baseline and recording should be done three times a day and at the same time extra measures should be taken if any changes in the health are observed



Result: The pulse oxygen (before and after exertion) was studied.

Hello

Friends

If you Get Any Help From This Notes /

Videos

Next You Turn To Help Us

Please Contribute Some Amount

To Our

FDSPharmacy Team

Phone pe 6398439940

Paytm 6398439940

Google Pay 6398439940



Amir Khan

Thank You Keep Supporting

